com-TEns 111

Transcutaneous Ele€trical Nerve Stimulation

Model: GM3AlOT/GM3AlOT(O)

Care for a Healthy Life

CE 0120

Please read the instruction manual befare use.

USER'S MANUAL

English

Introduction..... 1

Español

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Português

Introdução..... 42

SYMBOLS

Manufacturer

Authorised.representative.in.the.European.Community

CE.mark.and.notified.body.registration.number,.the.Annex.II.(excluding.section.4).article.3.of.EC.directive.93/42/EEC.have.been.met..

“BF”.symbol,.indicate.this.product.is.according.to.the.degree.of.protecting.against.electric.shock.for.type.BF.equipment..

GENERAL DESCRIPTION

Electric.myostimulation.has.proven.its.high.value.as.a.method.of.pain.therapy.and.is.a.great.help.to.the.experienced.therapist.

With.some.indications,.physicians.can.prescribe.unit.to.patients.for.the.use.at.home.

The.unit.is.a.dual-Channel.electric.stimulator.for.active.treatment.application,.which.output.as.well.as.an.8-bit.micro.computer.for.controlling.the.system.

The.electronics.of.the.unit.create.electric.impulses:.the.intensity,.duration,.frequency.per.second.and.modulation.of.these.impulses.can.be.adjusted.through.the.switch.or.knob.

SYSTEM COMPONENTS

Your.device.may.include.the.following.components.or.accessories:

- ***Unit***
- ***Carrying case***
- ***Lead wires / Electrodes***
- ***9-volt battery***
- ***Operation Manual***

WARRANTY

This.device.carries.a.one-year.warranty.from.the.date.of.purchase..The.warranty.applies.to.the.device.and.necessary.parts.and.labor.relating.thereto..The.distributor.reserves.the.right.to.replace.or.repair.the.unit.at.their.discretion.

The.warranty.does.not.apply.to.electrode,.battery,.lead.wires,.carrying.case,.damage.resulting.from.failure.to.follow.the.operating.instructions,.accidents,.abuse,.alterations.or.disassembly.by.unauthorized.individuals.

INDICATIONS AND CONTRAINDICATIONS

Read.the.operation.manual.before.using.the.device.

Federal.law.(USA).restricts.this.device.to.sale.by.or.on.the.order.of.a.physician..Observe.your.physician's.precise.instructions.and.let.him.show.you.where.to.apply.the.electrodes..For.a.successful.therapy,.the.correct.application.of.the.electrodes.is.an.important.factor.

Carefully.write.down.the.settings.your.physician.recommended.

Indications for use

This.device.is.a.prescription.device.and.only.for.symptomatic.relief.of.chronic.intractable.pain.

Contraindications

- . Any.electrode.placement.that.applies.current.to.the.carotid.sinus.(neck).region.
- . Patients.with.implanted.electronic.devices.(for.example,.a.pacemaker).or.metallic.implants.should. not.undertake.
- . Any.electrode.placement.that.causes.current.to.flow.transcerebrally..(through.the.head).The.use.of. unit.whenever.pain.symptoms.are.undiagnosed,.unit.etiology.is.determined.
- . The.use.of.TENS.whenever.pain.syndromes.are.undiagnosed,.until.etiology.is.established.

WARNINGS AND PRECAUTIONS

⚠Warnings

- The.device.must.be.kept.out.of.reach.of.children.
- The.safety.of.device.for.use.during.pregnancy.or.delivery.has.not.been.established.
- Do.not.place.electrodes.on.front.of.the.throat..This.may.result.is.spasms.of.the.laryngeal.and.pharyngeal.muscles.
- Do.not.place.the.electrodes.over.the.carotid.nerve.
- The.device.is.not.effective.for.pain.of.central.origin.(headaches).
- The.device.may.interfere.with.electronic.monitoring.equipment.(such.as.ECG.monitors.and.ECG.alarms).
- Electrodes.should.not.be.placed.over.the.eyes,.in.the.mouth,.or.internally.
- These.devices.have.no.curative.value.
- TENS.devices.should.be.used.only.under.the.continued.supervision.of.a.physician.
- TENS.is.a.symptomatic.treatment.and.as.such.suppresses.the.sensation.of.pain.which.would.otherwise.serve.as.a.protective.mechanism.

⚠Precautions/adverse Reactions

- Isolated. cases. of. skin. irritation. may. occur. at. the. site. of. electrode. placement. following. long-term. application.
- Stimulation. should. be. stopped. and. electrodes. removed. until. the. cause. of. the. irritation. can. be. determined.
- Effectiveness. is. highly. dependent. upon. patient. selection. by. a. person. qualified. in. the. management. of. pain. patients.
- If. the. device. treatment. becomes. ineffective. or. unpleasant,. stimulation. should. be. discontinued. until. reevaluation. by. a. physician/clinician.
- Always. turn. the. device. off. before. applying. or. removing. electrodes.
- Skin. irritation. and. electrode. burns. are. potential. adverse. reactions.


DANGER

The. device. does. not. have. AP/APG. protection.. Explosion. hazard. is. possible. if. used. in. the. presence. of. explosives,. flammable. materials. or. flammable. anesthetics.. Caution. should. be. used. when. applying. the. device. to. patients. suspected. of. having. heart. disease.. Further. clinical. data. is. needed. to. show. if. there. are. adverse. side. effects. on. those. with. heart. disease.


ABOUT THE DEVICE

Your. device. offers. two. controllable. output. channels.. This. device. creates. electrical. impulses. whose. amplitude,. width,. rate. and. modulation. can. be. altered. with. the. switches. or. knobs.. The. device. controls. are. very. easy. to. use. and. the. slide. cover. protects. accidental. changes. in. settings.

GM3A10T: With Timer

GM3A10T(0): Without Timer

THE DEVICE CONTROLS

Slide cover

A.cover.conceals.the.controls.for.Mode,.Time.(for.GM3A10T.only),.Width,.Rate..Press.the.topside.of.the.cover.and.pull.down.in.order.to.open.the.cover.

Intensity

The.intensity.knobs.are.located.on.the.top.of.the.unit.for.the.strength.adjustment.of.the.stimulation.and.also.function.as.ON/OFF.controls.

Mode

The.Mode.switch.is.used.to.select./.set.the.type.of.treatment.utilized..The.three.modes.are.Burst.(B),.Continuous.(C),.Modulation.(M).

TIME (for GM3A10T only)

Treatment.Time.of.device.can.be.pre-select./.set.with.Time.switch..There.are.two.programs.fixed.duration.of.30.and.60.minutes.and.one.program.of.continuous.output..Set.the.switch.to.the.position.desired.

To.reset.timer,.turn.both.channels.off...When.you.are.ready.to.begin.another.treatment.cycle.(at.least.rest.for.7.seconds),.turn.the.unit.on.as.normal...The.timer.will.be.reset.to.its.current.setting.

WIDTH

The.pulse.Width.controller.regulates.the.pulse.width.for.both.channels.

RATE

The.pulse.Rate.controller.regulates.the.number.of.pulse.per.second.for.both.channels.

ATTACHING THE LEAD WIRES

The.lead.wires.provided.with.the.device.insert.into.the.jack.sockets.located.on.top.of.the.unit..Holding. the.insulated.portion.of.the.connector,.and.pushing.the.plug.end.of.the.wire.into.one.of.the.jacks..After. connecting.the.wires.to.the.stimulator,.attach.each.wire.to.an.electrode.

Lead.wires.provided.with.the.device.are.compliant.with.mandatory.compliance.standards.as.set.forth. by.the.FDA.

Note: Use carefully when you plug and unplug the wires. Pulling on the lead wire instead of its insulated connector may cause wire breakage.

Caution: Never insert the plug of the lead wire into an AC power supply socket.

ELECTRODE SELECTION AND CARE

Your. physician/practitioner. should. decide. which. type. of. electrode. is. best. for. your. condition.. Follow. application.procedures.outlined.in.electrode.packaging.to.maintain.stimulation.and.prevent.skin.irrita- tion..The.electrode.packaging.will.provide.instruction.for.care,.maintenance.and.proper.storage.of.your. electrodes.

Be.sure.to.use.the.electrodes.provided.by.the.manufacturer.and/or.the.similar.FDA.legally.marketed. electrode.in.particular.the.same.cross.section.area.

TIPS FOR SKIN CARE

Good.skin.care.is.important.for.comfortable.use.of.your.device.

- Always.clean.the.electrode.site.with.mild.soap.and.water.solution,.rinse.well,.and.blot.dry.thoroughly. prior.to.any.electrode.application.
- Any.excess.hair.should.be.clipped,.not.shaved,.to.ensure.good.electrode.contact.with.the.skin.
- You.may.choose.to.use.a.skin.treatment.or.preparation.that.is.recommended.by.your.physician.. Apply,.let.dry,.and.apply.electrode.as.directed..This.will.both.reduce.the.chance.of.skin.irritation.and. extend.the.life.of.your.electrodes.
- Avoid.excessive.stretching.of.the.skin.when.applying.electrodes,.this.is.best.accomplished.by.applying.the.electrode.and.smoothly.pressing.it.in.place.from.the.center.outward.
- When.removing.electrodes,.always.remove.by.pulling.in.the.direction.of.hair.growth.
- It.may.be.helpful.to.rub.skin.lotion.on.electrode.placement.area.when.not.wearing.electrodes.

CONNECTING THE DEVICE

1. Prepare the Skin

Prepare.the.skin.as.previously.discussed.and.according.to.instructions.provided.with.your.electrodes.. Before.attaching.the.electrodes,.identify.the.area.in.which.your.physician/practitioner.has.recommended.for.electrode.placement.

2. Connect lead wires to the electrodes

Connect.the.lead.wires.to.the.electrodes.before.applying.the.electrodes.to.the.skin.

Note: Be sure both intensity controls for Channel 1 and 2 are turned to the "OFF" position.

3. Place Electrodes on Skin

Place.the.electrodes.on.the.skin.as.recommended.by.your.clinician.

4. Insert Lead Wire Connector to device

Plug.end.of.lead.wire.into.the.channel.output.receptacle.to.be.used,.pushing.plug.in.as.far.as.it.will.go.

5. Select Treatment Settings

Check.and.be.sure.your.unit.is.set.to.the.proper.settings.as.recommended.by.your.physician/practitioner.

6. Adjusting Channel Intensity Control

Locate.the.intensity.control.knobs.at.the.top.of.the.unit.under.the.cap..Open.the.cap.first,.slowly.turn.the.intensity.control.knob.for.Channel.1.clockwise.until.you.reach.the.intensity.recommended.by.your.medical.professional..Always.start.with.the.lowest.step.and.increase.slowly..Repeat.the.same.process.for.Channel.2,.if.appropriate..Close.the.cap.after.intensity.setting.ok.

If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation amplitude to a comfortable level. Cease stimulation and contact your physician/practitioner if problems persist.

BATTERY INFORMATION

A.9-volt.disposable.battery.is.provided.with.your.unit..Dimming.of.the.indicator.light.signifies.that.the.battery.should.be.replaced.with.a.new.one.as.soon.as.possible..However,.the.stimulator.will.continue.to.operate.for.several.more.hours.

CHANGING THE BATTERY

To.replace.the.battery:

1.. Remove.the.slide.cover.by.pressing.the.top.and.sliding.down.until.it.is.completely.removed.from.the.unit..This.will.reveal.the.battery.compartment.
2.. Remove.the.discharged.battery.from.the.device.
3.. Place.a.new.battery.in.the.compartment..Note.the.proper.polarity.alignment.indicated.on.both.the.battery.and.the.compartment.

CLEANING FOR YOUR DEVICE

Your.device.may.be.cleaned.by.wiping.gently.with.a.damp.cloth.moistened.with.mild.soap.and.water..Never.immerse.the.device.in.water.or.other.liquids.

Wipe.lead.wires.with.a.damp.cloth.as.above.if.they.become.soiled.

To.properly.store.the.device.for.an.extended.period.of.time,.remove.the.battery.from.the.unit..Put.the.unit.and.accessories.in.the.carrying.case.and.store.in.a.cool,.dry.location.

TROUBLESHOOTING

If.the.device.does.not.function.properly:

1.. Make. sure. the. battery. is. properly. installed,. or. replace. the. battery.. Be. sure. to. observe. proper. polarity. markings. when. replacing. the. battery.. If. the. indicator. light. is. not. on. when. the. unit. is. turned.on,.replace.the.battery.and.check.again.
2.. If. the. intensity. has. been. adjusted. and. there. is. no. stimulation,. check. that. the. lead. wires. are. properly.connected.and.the.electrodes.are.in.place..If.the.unit.appears.to.be.functioning.and.no. stimulation.occurs,.the.lead.wires.or.electrodes.may.need.to.be.replaced.
3.. If. the. battery. appears. to. be. charged. and. the. unit. is. not. functioning,. turn. both. intensity. Control. Knobs.to.the.OFF.position.(counter-clockwise).for.about.7.sec..Then.gradually.turn.the.intensity. Control.Knob.clockwise.until.stimulation.is.felt..If.device.still.is.not.working,.turn.the.unit.off.and. contact.your.distributor.

If.any.other.problems.occur,.please.consult.or.return.the.device.to.your.distributor..Don't.try.to.repair. a.defective.device.

TECHNICAL SPECIFICATIONS

Channel	**Dual, isolated between channels**
Pulse intensity	**Adjustable 0-80mA peak into 500 Ω load each channel, constant current**
Pulse Rate	**2Hz-120Hz (adjustable).**
Pulse Width	**40μs - 260μs (adjustable).**
Timer (for GM3A10T only)	**30, 60 minute and continuous mode selectable**

Function Modes:

TENS	
B	**Cycle Bursts, 2Bursts/sec, 9pulses/Burst, 100Hz, width is adjustable.**
C	**Continuous mode. Pulse rate, pulse width and intensity are adjustable.**
M	**Modulated Width. Pulse width is automatically varied in an interval of 6 seconds.The modulation range of pulse width is from setting value to 35% less than the control setting value, then returns to the setting value. Rate, width and intensity are fully adjustable.**
Wave form	**Asymmetrical Bi-phasic square pulse.**
Voltage	**0-110 Volt (Open Circuit).**
Max charge per pulse	**21 micro-coulombs**
Power Source	**9-Volt Battery.**
Dimensions	**95 mm(H) ×61.5mm(W) ×26m(T).**
Weight	**Approx. 120grams (battery included).**

All electrical specifications are ±10% except the amplitude is +/- 20% (500Ω load).

EMC INFORMATION			
Recommended. separation. distances. between. portable. and. mobile. RF. communications. equipment. and. the. ME. equipment			
The.ME.equipment.is.intended.for.use.in.an.electromagnetic.in.which.radiated.RF.disturbances.are.controlled..The. customer.or.the.user.of.the.ME.equipment.can.help.prevent.electromagnetic.interence.by.maintaining.a.minimum. distance. between. portable. and. mobile. RF. communications. equipment. (transmitters). and. the. ME. equipment. as. recommended.below,.according.to.the.maximum.output.power.of.the.communications.equipment.			
Rated.maximum.output.power.of.transmitter W	Separation.distance.according.to.frequency.of.transmitter.m		
	150 kHz to 80 MHz d=1,2 $\sqrt{P}$	***80 MHz to 800 MHz*** d=1,2. $\sqrt{P}$	***800 MHz to 1.0 GHz*** d=2,3 $\sqrt{P}$.
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

Declaration – electromagnetic emissions and immunity –

for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

The ME declaration – electromagnetic immunity			
The.ME.system.is.intended.for.use.in.the.electromagnetic.environment.specified.below.. The.customer.or.the.user.of.the.ME.system.should.assure.that.it.is.used.in.such.an.environment.			
Immunity.test	IEC. 60601. test. level	Compliance. level	Electromagnetic.environment.-.guidance.
Conducted.RF IEC.61000-4-6 Radiated.RF IEC.61000-4-3	3.Vrms 150.kHz.to.80.MHz 3.V/m 80.MHz.to.1.GHz	3V 3V/m.	Portable. and. mobile. RF. communications. equipment. should.be.used.no.closer.to.any.part.of.the.EQUIPMENT. or. SYSTEM. including. cables,. than. the. recommended. separation. distance. calculated. from. the. equation. appli-cable.to.the.frequency.of.the.transmitter. Interference.may.occur.in.the.vicinity.of. equipment.marked.with.the.following.symbol. 
Declaration – electromagnetic immunity			
The.ME.system.is.intended.for.use.in.the.electromagnetic.environment.specified.below. The.customer.or.the.user.of.the.ME.system.should.assure.that.it.is.used.in.such.an.environment.			
Immunity.test	IEC. 60601. test. level	Compliance. level	Electromagnetic.environment.-.guidance
Electrostatic.dis-charge.(ESD) IEC.61000-4-2	±3.kV.contact ±8.kV.air	±3.kV.contact ±8.kV.air	Floors.should.be.wood,.concrete.or.ceramic.tile..If.floors.are. covered.with.synthetic.material,.the.relative.humidity.should. be.at.least.30.%.

Declaration – electromagnetic emissions		
The.ME.system.is.intended.for.use.in.the.electromagnetic.environment.specified.below. The.customer.or.the.user.of.ME.should.assure.that.it.is.used.in.such.an.environment.		
Emissions.test	Compliance	Electromagnetic.environment.-.guidance

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Print-2012/All rights reserved